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de

Pe-Ha-Visco®

Hydroxypropylmethylcellulose

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Hydroxypropyl methylcellulose

1. Description

The product is a clear, sterile, isotonic and viscoelastic preparation of hydroxypropyl methylcellulose (HPMC) for intraocular and extraocular use, dissolved in a buffered physiological solution (pH 6.8-7.5). The hydrogel is filled in a ready-to-use syringe for single use, which is packaged in a sterile blister pack. Depending on the product version it contains the following concentration of HPMC: 2.0% or 2.4%. The product is free from preservatives. It has been subjected to moist heat sterilization for injectable products and is non-pyrogenic. There are no known inflammatory or immunogenic reactions. Chemical composition (w/v) as mg of substance per 1 ml hydrogel: See table.

2. Product characteristics and mode of action

As an aid during surgery of the anterior segment of the eye, the product maintains the depth of the anterior chamber during intraocular surgery. The product supports intraocular surgery thanks to its viscoelastic and wetting properties. Due to its water solubility, the product is easily removed by irrigation and aspiration at the end of surgery.

3. Indications

The product is used as a volume substitute for the aqueous humor during intraocular surgery such as lens extraction and insertion of an intraocular lens. The product maintains the depth of the anterior chamber throughout surgery and reduces the risk of traumatization of the corneal endothelium, the iris and the ciliary body as a result of direct contact with surgical instruments. When applied extraocularly, the product acts as a wetting agent during ocular surgery and as a lubricant and wetting agent during the implantation of intraocular lenses using injector or cartridge systems.

4. Contraindications

The product should not be used in patients with known hypersensitivity to hydroxypropyl methylcellulose or to any other constituent of the hydrogel. Children, pregnant or nursing women should not be treated with hydroxypropyl methylcellulose as there is no clinical data available on its use in these patients.

5. Precautions

Observe all routine precautions which are required during ophthalmic surgery. The product, together with the enclosed single-use cannula, sterilized using ethylene oxide, is intended for single intraocular or extraocular use in one eye only. It must not be re-sterilized, resealed or reused. Do not re-sterilize the solution, as this can cause changes to the product characteristics. Do not use the product if there is any damage to the ready-to-use syringe or sterile packaging. The product is to be used for one patient only during a single appointment.

6. Potential side effects

Biocompatibility testing of the product did not show undesirable effects. After the intraoperative use of the product, a transient post-operative rise in intraocular pressure may occur. Suitable treatment to reduce the intraocular pressure can counteract this increase.

The healthcare professional should advise the patient

- of any potential adverse events related to the product.

- that the patient must report any adverse event or complication to a physician.

Serious events that have occurred in connection with the use of the product must be reported to the manufacturer and to the competent national authority in which the physician and/or patient is established.

7. Drug and chemical interactions

To date, there is no available data on the incompatibility of the product with other products for use on or in the eye.

8. Dosage, mode and duration of treatment

The enclosed sterile ophthalmological single-use cannula (or a comparable sterile cannula) is to be screwed onto the luer lock connector of the syringe and vented before the product is injected into the eye. For intraocular lens implantation, coat the implant and the instruments with the product immediately before surgery. This also protects the endothelium and surrounding tissue. The volume of product for injection will vary based on the individual patient and the surgical procedure. To compensate for loss of viscoelasticity caused by leakage or irrigation, several injections of the product may be performed. Do not exert any further pressure on the plunger prior to removing the cannula from the eye. This prevents the aspiration of air bubbles into the cannula. Remove all product at the end of surgery using a suitable irrigation/aspiration device.

9. Dosage form

a) Box with one sterile ready-to-use syringe containing 2.0 ml hydrogel and one sterile single-use luer lock cannula.
b) Box with 10 sterile ready-to-use syringes containing 2.0 ml hydrogel and 10 sterile single-use luer lock cannulas.

10. Shelf life

Do not use the product after the expiry date. The expiry date (year/month) is printed on the syringe blister pack and box.

11. Storage

Store in the original packaging in a dry place at room temperature and away from light. Avoid freezing or impact. Always observe the information symbols on the packaging. Store this product in a location inaccessible to children.

12. User information

If a product from damaged packaging is used, or a product reused on other patients, or on the same patient at a later time, the intended properties and sterility of the product cannot be guaranteed. Potential material contamination and transfer of pathogens from the general environment, the clinical environment or (other) patients could lead to health complications. Resealing or re-sterilization is not a permitted or guaranteed method of restoring the product to a usable condition.

13. Recommended user type

For use solely by trained ophthalmologists.

14. Disposal

The syringe contents are non-toxic and non-flammable. Unused syringes and their contents are not infectious and may be disposed of before or after their expiry date, in accordance with national and local regulations. Used syringes and cannulas must be disposed of as epidemiologically hazardous waste in line with national and local regulations for safe handling and disposal.

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