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Package leaflet and prescribing information **Please read carefully!**

Tutoplast® Fascia lata Tissue Matrix

- Name of the medicinal product**
Tutoplast® Fascia lata
- Qualitative und quantitative composition**
Human fascia lata, preserved using Tutoplast® technique, sterilised with gamma irradiation.
The package contains one graft.
- Pharmaceutical form**
Tissue graft
- Clinical particulars**
4.1 Therapeutic indications
As a connective tissue substitute or graft for tension-free or unidirectional loading as a tendon, ligament or fascia replacement among other uses. Positive experiences have also been reported for the clinical use of the product in neurosurgery, orthopaedics, trauma surgery, maxillofacial surgery, ophthalmology, urology, gynaecology, plastic surgery and ENT surgery.

4.2 Posology and method of administration
The dimensions of the graft are determined by the operating surgeon based on the defect to be repaired. It may prove necessary to cut the membrane to ensure a per-

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fect fit. The graft must be completely rehydrated in normal saline solution or Ringer's solution prior to use (see 6.6). The transplantation for permanent integration in the human body is performed once the graft is completely rehydrated.

- Contraindications**
None known
- Special warnings and precautions for use**
 - Strict assessment of the indications must be performed in the following cases prior to use of the graft:
 - For implantation in necrotic host sites
 - For implantation in regions with poor blood supply
 - For implantation in host sites with an active or latent infection
 - In the case of disorders or conditions which could have a negative impact on the healing rate
 - Based on their medical history, as well as further tests, all tissue donors are examined for possible grounds for exclusion. Before the donor tissue can be approved for processing, serological tests are carried out for possible infectious agents. This analysis includes tests for HTLV-I/-II, HIV-1/-2, hepatitis B/C and Treponema pallidum, among others. The Tutoplast® process is effective against all known pathogenic species. However, as with all biological products, the transmission of infectious diseases can be ruled out only to the extent of current knowledge
 - The Tutoplast® graft remains sterile as long as the packaging is undamaged. If the sterile packaging becomes damaged, the product must be disposed of appropriately. Do not resterilise the product! If the graft is contaminated in the course of a surgical intervention, it must be disposed of appropriately.
 - Correct positioning and securing in place of the graft (see also 6.6.) are essential prerequisites for its successful application. The graft must be implanted and secured in such a way that it remains in the intended position and its incorporation is thus facilitated.

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- Tutoplast® Fascia lata should only be used to cover defects which are completely free from tension or only subjected to unidirectional stress. The tissue can withstand high loads along the grain of the fibre. However, mechanical loading across the grain of the fibre is contraindicated. It is recommended that the fascia be doubled at suture points subject to particular stress.
 - Tutoplast® Fascia lata is intended for single use only. Unused material must be disposed of appropriately.
- 4.5 Interaction with other medicinal products and other forms of interaction**
None known
- 4.6 Pregnancy and lactation**
No clinical data concerning the use of the product during pregnancy and when breastfeeding are available. Based on the physical mechanism of action of the graft, undesirable effects during pregnancy or on the health of the newborn child are not to be expected and have not been determined.
- 4.7 Effects on ability to drive and use machines**
Not relevant
- 4.8 Undesirable effects**
The following classifications are used as the basis for the frequency of adverse reactions:
- Very common (≥1/10)
Common (≥1/100 to <1/10)
Uncommon (≥1/1,000 to <1/100)
Rare (≥1/10,000 to <1/1,000)
Very rare (<1/10,000)
Not known (frequency cannot be estimated from the available data)

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System organ class	Frequency	Adverse effect
Immune system disorders	Not known	Graft rejection
General disorders and administration site conditions	Not known	Reaction at the site of implantation
Injury, poisoning and procedural complications	Not known	Graft failure

As with any operation/graft, there is also a risk of infection or another reaction caused by the intervention.

- Reporting of suspected adverse reactions**
Reporting suspected adverse reactions of the medicinal product after its approval is important.
It allows continued monitoring of the benefit/risk balance of the medicinal product.
• Healthcare professionals within Germany are asked to report any suspected adverse reactions to the Bundesinstitut für Impfstoffe und biomedizinische Arzneimittel (Federal Institute for Vaccines and Biomedicines), Paul-Ehrlich-Institut, Paul-Ehrlich-Str. 51-59, 63225 Langen, telephone: +49 6 10 37 70, Telefax: +49 61 03 77 12 34, email: biovigilance@pei.de, website: www.pei.de.
• Healthcare professionals outside Germany are asked to report any suspected adverse reactions to the manufacturer: Tutogen Medical GmbH, Industriestrasse 6, 91077 Neunkirchen am Brand, Germany Tel.: +49 9134-9988-0, Fax: +49 9134-9988-99, vigilance@rtix.com National legislation may stipulate additional requirements for reporting to local authorities.

Patients must be informed to contact their doctor or medical professional if they notice any side effects. This also applies to side effects that are not stated in this leaflet and summary of product characteristics. Patients can also report side effects

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directly to the Paul Ehrlich Institute. By reporting side effects, patients can help provide more information about the safety of this medicinal product.

4.9 Overdose

Not relevant

5. Pharmacological properties

5.1 Pharmacodynamic properties

Pharmacotherapeutic group: Allograft

Mechanism of action of the medicinal product

Tutoplast® Fascia lata ensures reliable closure of the tissue defect and serves as a guide rail for the cellular repair mechanism. The regeneration of autologous structures occurs in parallel to the resorption of the graft. In addition to on the size of the graft itself, both processes are also dependent on the blood supply to the surrounding tissue. The ultimate goal is complete conversion of the graft into autologous tissue.

5.2 Pharmacokinetic properties

Tutoplast® Fascia lata is primarily a physical aid, acting as a guide rail for the regeneration of soft tissue. The process of cell colonisation, resorption and remodelling generally begins after 1-2 days. The conversion process depends on the size of the graft and the host site's ability to respond. As such, it can take a considerable period of time to complete (from several weeks to several years).

5.3 Preclinical safety data

Studies on animal models and *in vitro* studies have shown that the graft is biocompatible, is incorporated quickly and maintains the space until replaced by autologous tissue. As such, the preclinical data do not indicate any risks to humans and confirm that the product can be used safely.

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None known

6.3 Shelf life

Tutoplast® Fascia lata has a shelf life of 5 years providing the original packaging is intact. The graft must not be used after the expiry date shown on the packaging.

6.4 Special precautions for storage

The graft should be stored in a clean, dry place protected from exposure to direct sunlight and at no higher than 30°C. The product must not be frozen.

6.5 Nature and contents of container

The graft is kept sterile in a double-layered primary packaging. The packaging contains a pre-cut piece of Tutoplast® Fascia lata.

	Shape, size (mm)	Packaging unit
Tutoplast® Fascia lata Tissue Matrix	3 x 120, 3 x 140, 10 x 30, 10 x 100, 10 x 140, 15 x 50, 15 x 60, 15 x 200, 15 x 300, 20 x 30, 20 x 40, 20 x 70, 20 x 100, 20 x 120, 20 x 140, 20 x 190, 20 x 200, 20 x 240, 20 x 300, 25 x 30, 30 x 40, 30 x 60, 40 x 50, 40 x 70, 40 x 100, 40 x 200, 50 x 110, 60 x 80, 60 x 120, 70 x 20, 70 x 80, 70 x 120, 80 x 200, 200 x 40, Ø 14, 15 x 60, thickness ≤ 0.6, 15 x 60, thickness > 0.6	1 piece

Note: It is possible that not all packaging sizes will be available on the market.

6.6 Special precautions for disposal and other handling

Information on disposal of remnants

If the packaging is opened and the graft not used, the fascia lata must not be resteri-

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lised. Remnants must be disposed of in accordance with the pertinent national regulations. The contents of the packaging are exclusively intended for use in **one** patient.

Special information on handling

- Please observe the general principles of aseptic handling when using Tutoplast® Fascia lata.
- Before use, the graft must be completely rehydrated in normal saline solution or Ringer's solution (until the membrane is soft and flexible and can be well moulded to the surface).
- In the case of rectangular fascia pieces, the longer side always follows the grain. The size of the implant depends on the size of the respective tissue defect.
- The implantation should be performed in such a way that the free edges of the fascia do not extend into the abdominal cavity, for example (risk of adhesion).
- Autologous blood, blood components or bone marrow aspirate may only be added after rehydration and immediately prior to the product's implantation.
- Tutoplast® Fascia lata can be fixed in place with absorbable or non-absorbable sutures.
- Suture material with a round needle in an atraumatic version is selected for the continuous suture. The suture thickness depends on the indication and operation.
- The sutures should be placed on the graft at least 2-3 mm from the edge.

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Each graft is supplied with stickers for documentation in the patient's file. Tutoplast® is a registered trademark of Tutogen Medical GmbH.